ADVARSLER

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WARNINGS

Read the manual in entirety before using the apparatus.

All work with high voltages is potentially highly dangerous and must only be performed by persons qualified.

If you choose to connect this equipment to a high voltage, it is your responsibility to work in accordance with applicable rules and in compliance with all necessary precautions.

Do not connect the two plugs to voltage differences larger than 10 kV AC.

Make sure that the high voltage supply “floats” relative to ground.


Sparks and electric arcs will form ozone and nitrous gases.

Ensure ventilation

Turn off the high voltage supply (at the mains) before touching the plugs. This applies to both connection and disconnection.

Touch the plugs with one hand only (keep the other on the back.)

Never turn on the high voltage unless **both** plugs are inserted.

Avoid overheating of the high voltage transformer used.
Don't let the experiment run for more than 2 to 3 minutes.

The distance between the electrodes can be varied using the black plastic rod at the bottom of the device.

When voltage is applied, do not touch any other part of the apparatus.

Touch the outermost part of the rod with one hand only (keep the other on the back.)

In case of any sign of damage to the device, including cables and plugs, do not use it.

Formation of an electric arc

When the electrodes are close enough to each other, a high voltage will cause a spark to jump between them. The spark ionizes the air molecules, resulting in a sustained discharge.

The electrical energy deposited will heat the air, causing it to move upwards. At some point, the length of the arc – and thus the electrical resistance – will become so large that the discharge stops.

High voltage supply

A transformer with a primary coil with 600 windings and a secondary coil with 24000 windings can transform a primary voltage of 230 V to a secondary voltage of 9.2 kV.

Frederiksen Scientific can provide the following items:

- 459645 Coil 600 wdg. For mains (SCHUKO plug)
- 459655 Coil 24000 wdg. Safety sockets
- 459700 UI-core with clamping brackets

(Laboratory high-voltage supplies – like for instance Frederiksen's 367060 – are usually provided with a current limiter which puts a very low limit on the output power.

Such a power supply is completely unsuitable for the production of powerful electric arcs.)

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© Frederiksen Scientific A/S