2014.11.28

4565.10 AC

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4565.10 AC

## Description of the equipment:

This equipment is designed to measure the force acting upon an electric conductor in a magnetic field, compared to the angle between these. When an electric conductor with the length  $L$  is within a magnetic field  $B$ , and a current  $I$  is running through the conductor, this will be acted upon by a force, which can be calculated by the following formula (also known as Laplace's law):

$$F = B \cdot I \cdot L \cdot \sin(V)$$

Here  $V$  is the angle between the current and the magnetic field.

The equipment consists of a holder with permanent magnets, ensuring an almost homogenous magnetic field, and a revolving coil in a goniometer assembly, where one can read the angle between the current and the magnetic field. The goniometer assembly fits into the holder for wireframes from current balance 1.

With this equipment the relation between force and angle can be determined. The force is measured indirectly. The magnet holder is placed on a balance, which you tare. The force is found by applying Newton's third law. In this situation it means that the

force acting upon the conductor from the magnetic field is equal in magnitude, but opposite in direction, to the force acting upon the magnetic field from the current in the conductor. If the force acting on the conductor is directed upwards, then the force acting on the magnetic holder will be directed downwards, which results in an apparent increase of the weight of the magnet holder. Since the force is proportional to the weight ( $F = m \cdot g$ ), the force can easily be calculated.

For a more detailed description of the experiment, please refer to the experimental section of current balance 1.

This can be downloaded from [www.frederiksen.eu](http://www.frederiksen.eu)

## Necessary accessories:

Current balance 1, type 4565.00

Powersupply with stabilized DC, 3630.00 or 3640.00

Laboratory Balance with 10 mg resolution and at least 200 g capacity.

Base tripod 0006.00

Rod 0008.50

Test cables

## Operation:

The goniometer assembly is mounted on the holder for wireframes, which comes from current balance 1. Be careful to adjust the goniometer so that the moveable angle indicator reads zero, when the horizontal (bottommost) conductors of the coil are parallel to the magnetic field. Tare the sensitive balance with no current flowing. Now supply a constant current of e.g. 4 amperes, and adjust the angle between the conductors and the magnetic field in 10 degree increments. Read off corresponding "force" values from the balance.

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© A/S Søren Frederiksen, Ølgod